

Sheet2

<i>values as found</i>	<i>status as found</i>	<i>changed?</i>	<i>status, current</i>
isct1 als cavity	works	no change	works
SYS-MOTION_C_SHUTTER_D_NAME			
isctex green beam	works	no change	works
SYS-MOTION_X_SHUTTER_A_NAME			
isctex fiber beam	appears to not work	no change	appears to not work
SYS-MOTION_X_SHUTTER_B_NAME			
isct6 spare	appears to not work	no change	appears to not work
SYS-MOTION_C_SHUTTER_H_NAME			
isctey green beam	appears to not work	changed to	works
SYS-MOTION_Y_SHUTTER_A_NAME		SYS-MOTION_Y_SHUTTER_B_NAME	
isctey fiber beam	works, but shuts off the isctey green beam	changed to	appears to not work
SYS-MOTION_Y_SHUTTER_B_NAME		SYS-MOTION_Y_SHUTTER_A_NAME	