
















DIAGNOSTIC PD

DC OUT

NE GUARD

MOPA DIAG

INVISIBLE LASER RADIATION
AVOID DIRECT EXPOSURE TO BEAM
MAXIMUM OUTPUT POWER 3000mW
WAVELENGTH 1064 nm
USE EYE PROTECTION
CLASS IV LASER PRODUCT

This product complies with 21 CFR 1040.10

INNOLIGHT 

Model:

Mephisto 2000 NE LIGO


Serial No: 1639E

Made in Germany

Manufactured: 10 / 2006

INNOLIGHT GmbH, Garbsener Landstr. 10
D-30419 Hannover, Germany
www.innolight.de

Enhanced LIGO
Power Amplifier
PA-1208-7



This product complies with
21 CFR 1040.10 and
1040.11

Enhanced LIGO
MOPA System

Spare 1-1208-7

NOTICE

**NPRO 1639E inside
but not installed.**

Crystal temperature: 25.00 degC
Diode A temp: 31.62 degC, Diode B temp: 32.70 degC
injection current: 2.30 A
Output power: 2.09 W
PZT tuning: 1.1 MHz/V

NOTICE

NPRO 1639E inside
but not installed.

Crystal temperature: 25.00 degC
Diode A temp: 31.62 degC, Diode B temp: 32.70 degC
Injection current: 2.30 A
Output power: 2.09 W
PZT tuning: 1.1 MHz/V

DIAGNOSTIC PD

DC OUT

NE GUARD

WOPA DIAG
SERVO-POWER

COM 2

COM 1

ADW

POWER IN +/- 24

This product complies with
FCC CFR 47.101.10 and
47.101.11

External L&S
NPRO System
Serial 1.000-7

Serial Number 888220